

Zeolite MED[®] skin powder Instructions for use

ZEOLITE MED[®] SKIN POWDER FOR LOCAL APPLICATION TO THE SKIN ACCELERATES HAEMOSTASIS, SUPPORTS THE HEALING OF WOUNDS

Das Zeolite MED[®] skin powder is a non-invasive medical product and a natural bio-regulator for external use on the skin.

The medical product Zeolite MED[®] is made up of the natural mineral clinoptilolite zeolite; the proportion can reach up to 100%. This is a crystalline micro porous mineral of volcanic origin. Thanks to its large inner surface area and its extraordinary properties for selective ion exchange and absorption, Zeolite MED[®] can, where required, absorb and excrete exudates, harmful substances, in particular the heavy metals mercury, lead and cadmium, radionuclides, histamine and ammonium, in its micro porous structure.

These properties mean that Zeolite MED[®] skin powder can be applied, e.g. to damaged skin, as a protective shield for the wound; it can accelerate haemostasis, reduce the rate of microorganisms entering the wound and thus support the healing process.

AREAS OF APPLICATION AS A WOUND, BODY AND FOOT POWDER

- For detoxification of damaged skin, particularly grazes, light burns, insect bites, skin irritations
- For accelerated haemostasis
- For prevention of infection and rapid healing of wounds

WHEN SHOULD ZEOLITE MED[®] SKIN POWDER NOT BE APPLIED?

The skin powder should not be used on deep wounds that penetrate the dermis, which cannot be healed using only a secondary healing method. It is advisable to always seek medical care concerning deep, open wounds. If in doubt, please consult a doctor before application.

DOSAGE AND DURATION

Unless otherwise prescribed by your therapist, Zeolite MED[®] skin powder can be applied several times a day to the relevant body parts as a wound, body and foot powder, or to the area around smaller wounds. If necessary, the wound can be covered with a sterile dressing.

INTERACTION WITH OTHER SUBSTANCES

A waiting period of at least 2 hours should be observed following the use of medical ointments, creams etc.

WARNINGS AND CAUTIONARY MEASURES

The medical product is not a substitution for an antiallergic and cannot replace medical care of deep, open wounds! If in doubt, please consult a doctor before application.

The medical product must not be inhaled (breathed in) or used on the eyes. Should it come into contact with the eyes, rinse thoroughly.

Terminate the application if you notice an abnormal effect.

STORAGE AND USABILITY

Keep the product out of the reach of children aged 6 years and below.

Always close the package after use.

Store the product in a dry place at room temperature

The expiry date is printed on the bottom of the package. Do not use the product after the expiry date.

ADDITIONAL APPLICATION POSSIBILITIES

- Mix the package contents with water to form a paste, apply to the affected areas, cover with a compress and leave, e.g. overnight
- For use as a skincare product, take a pinch of Zeolite MED[®] skin powder in your hand, mix it with skin cream and apply as normal
- For masks, mix 1 tsp of skin powder with a small amount of water, apply in a thin layer, then rinse after 10 to 15 minutes

CONTENTS

Clinoptilolite zeolite

INSTRUCTIONS FOR DISPOSAL

Content and packaging can be disposed of with household waste.

MEDICAL PRODUCT MANUFACTURER

ZeoBent Handels GmbH

Weißer 2

D-07407 Uhlstädt-Kirchhasel

Germany

Phone +49 (0) 36742 673937

e-Mail: info@zeolith-bentonit-versand.de

Web: www.zeolith-bentonit-versand.de

<https://www.zeolite-bentonite.co.uk/en/zeolite-skin-powder-40g.html>



Status: 04/2018